



Days Bed Assist Handle

Item# 0816-10-815

Helping People Maintain Their Independence!



Soft Grip Handle
assists in rolling, sitting and getting out of bed

Pocket Organizer
helps keep items within arms reach

No Tools Required
quick and easy to put together on most beds

Steel Powder Coated Frame
supports up to 250 lbs

Height Adjustable Legs
handle height adjusts from 31" to 44"

Full Length Legs and Straps
provides added safety, stability and improves in bed mobility



Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200
©2013 Patterson Medical Products, Inc.

Fax toll free:
1-800-547-4333
In Canada fax 1-905-566-9207
081615624 LIT

Visit us online at:
www.pattersonmedical.com

