




TREATMENT/MAT TABLE, BARIATRIC



Motorized Hi-Lo, Scissor
Raised Back & Leg

<u>Item No.</u>	<u>Width</u>	<u>Length</u>
TWG8332.BTM	32"	78"
TWG8336.BTM	36"	78"
TWG8342.BTM	42"	78"
TWG8348.BTM	48"	78"

FEATURES:

- 750 lb Weight Capacity, *lifting*
- Motorized Height Adjustable
- Raised Back and Leg Sections, Motorized
- 4 Total Lock Casters
- 0 Horizontal Travel Space
- Hand Control
- Steel Top with 2" Removable Mat
- Steel Frame

Height Range: 17 1/2" - 36 3/4"
Finish: Pearl Grey Standard; Porcelain Optional.
Electrical: 120V, 6.5A, 780W, hospital grade plug
Warranty: Ten (10) Year Limited Warranty.

*Nominal Dimensions Stated