

# ROLYAN STABILIZER WALKER

## Product Features:

- Plastic molded uprights with spring steel reinforcement for increased durability.
- Contoured strut design allows the walker to conform to the leg, ensuring a better fit.
- Shock absorbing insole reduces impact of heel strike when ambulating.
- Low profile, rocker bottom is designed to promote a natural gait.
- Wide foot bed provides for increased comfort and stability.
- Lightweight construction.

## Indications:

- Acute ankle sprains, soft tissue injuries, stable fractures, post-operative use.



Low



High

Product #	Size	Men's Shoe Size	Women's Shoe Size
5500-11-01	XSmall	2 - 4	3 - 5
A414-7	Small	4.5 - 6	5.5 - 7
A414-8	Medium	6.5 - 10	7.5 - 11
A414-9	Large	10.5 - 12	11.5 - 13
5500-11-02	XLarge	12 +	13+

Product #	Size	Men's Shoe Size	Women's Shoe Size
9293-86-01	XSmall	2 - 4	3 - 5
A414-4	Small	4.5 - 6	5.5 - 7
A414-5	Medium	6.5 - 10	7.5 - 11
A414-6	Large	10.5 - 12	11.5 - 13
9293-88-01	XLarge	12 +	13+