

Euro Deluxe Chair™

0816-81-089 Shower Commode Chair



• Cushion Push/Pull Handle

• Tight Weave Mesh Back Support
Is Extremely Comfortable

• Cushion Infill Storage Pouch

• Foldaway Arms/
Removable Arms

• Lap/Chest Belts

• Soft Seat Cushion

• Commode Closure Cushion Infill

• Aluminium High Strength
Frame/Rust Proof

• Total Locking 4" Casters
2 Year Warranty

• Calf Straps

• Adjustable Height Foldaway
Footrest/Removable Footrest

• Cushion Infill Storage Pouch
(Storage Pouch will Help
Prevent from Mis-Placing)



• Durable and User-Friendly
Commode Bucket Comes with
Handle and Full Lid. Easy to Use,
Commode Pail is Held in Position
with Rails.

• Adjustable Height Frame

• Stainless Steel Fittings and Fixings
• 2 Year Limited Warranty on Frame Only

Euro Chair Specifications

• Weight Capacity	400lbs
• Seat Width	22"
• Seat Depth	17"
• Footrest Length	16" - 20"
• Backrest Fabric Height	13.7"
• Backrest Handle Height	19"
• Product (Net) Weight	30lb

• Arm Height from Seat	8.7"
• Dimension Between Arms	22"
• Dimension Between Vertical Legs (Internal Frame)	18.5"
• External Width	25.8"
• External Length (Footrests Included)	35.4"
• Overall Length (Footrests Excluded)	27.5"

• Seat to Floor Heights	1	19.6"
	2	20.6"
	3	21.6"
	4	22.6"
Clearance Heights	1	17.1"
	2	18.1"
	3	19.1"
	4	20.1"


PATTERSON
MEDICAL

www.pattersonmedical.com

Call toll free: **1-800-323-5547**
In Canada Call 1-800-665-9200

Fax toll free: **1-800-547-4333**
In Canada Call 1-877-566-8611