

# circaid® juxtait® essentials lower extremity

Inelastic, ready-to-wear compression system for the entire leg.

**Indications:**


- Mild to severe lymphedema.
- Mild to severe lipedema.
- Moderate to severe venous disease.

**Features:**

- Inelastic compression.
- Multiple compression levels (20-30, 30-40, 40-50 and 50+ mmHg).
- Juxtalock™ band system allows for easy and instant adjustability of product without removal.
- circaid® Built-In-Pressure system™ (BIPS™) guide card provides targeted compression ranges for improved results.
- SILV Eriect™ added to prevent static, odor & propagation of bacteria in the garment.
- Active & sedentary use.
- Simple measuring & sizing – Ready made in multiple sizes with 2 lengths of the lower leg and 3 lengths of the upper leg with knee.
- Works great with elastic stockings as part of a complete compression system.
- 6-month warranty.
- Available in beige.

**System includes:**

- 1 legging, 1 pair of knee high or high high under-sleeves (footless), 1 black cover up and 1 Built-In-Pressure system™ (BIPS™) guide card.
- Lower leg also includes: 1 circaid® power added compression band™ (pac band™) and 1 pair of compression anklets.
- See pages 40-42 for foot options.




juxtait® essentials upper leg with knee




juxtait® essentials lower leg with compression anklelet

sizing




circaid® juxtait® essentials lower leg											
circumference in centimeters											
size	S	MF	M-X full calf	L	L-X full calf	XL	XL'	XL-X full calf	XXL		
cC	30-38	36-43	43-50	41-48	48-56	41-51	46-56	51-61	53-64		
cB1	20-28	25-33	33-40	30-38	38-46	33-41	36-46	41-51	43-53		
cB	15-23	20-28	20-28	25-33	25-33	30-38	30-38	30-38	33-41		
actual product length in centimeters											
size	short					long					
l B-D	28					36					
sizing measurement											
l mall-D1	short					long					
	<36					>36					

\*Applies to circaid® readyfit lower leg and circaid® classicflex lower leg with adjustable lengths. Product lengths may differ from body measurements.



circaid® juxtait® essentials upper leg with knee												
circumference in centimeters												
size	XXS	XS	S	M	L	XL						
cG	40-52	46-62	54-70	60-80	69-90	77-101						
cE1	32-42	37-48	42-55	50-65	58-76	66-87						
cD	23-30	27-35	31-41	35-46	40-52	44-58						
actual product length in centimeters												
size	x-short			short			long					
l D-G	35			45			55					
sizing measurement												
l E-G	x-short			short			long					
	<25			<35			>35					

Product lengths may differ from body measurements.



circaid® power added compression band™ (pac band™)				
circumferences in cm.				
	standard	large	xlarge	
cB	19-27	25-50	<56	
cA1	19-36	22-41	<44	

\*circaid® compression anklelet only available in standard and large.

Leg Measurement Guide	
T	Waist
H	Widest part of pelvis
K	Gluteal fold
G	Groin
F	Mid-thigh
E2	15cm above E
E1	5cm above E
E	Center of patella
D1	Knee crease
D	Slightly below knee
C	Widest part of calf
B1	Between ankle and widest part of calf
B	Narrowest part of ankle above malleolus
A1	Middle of foot
A	Circumference of ball of foot
Y	Diagonally around heel over widest part of top of ankle
Z	Heel to base of great toe
mall	Center of malleolus
l	l-length
c	c-circumference