



VASYLI[®] MEDICAL & VIONIC[®]

ENHANCED
BIOMECHANICAL
FUNCTION



HELP RESTORE NATURAL FOOT FUNCTION
WITH **VASYLI ORTHOTICS & VIONIC FOOTWEAR**

Patterson Medical
EXCLUSIVE


VIONIC[®]

Natural Relief.

Footwear & orthotics help restore natural foot function and well-being with every step.

Our feet were designed to walk on soft, natural surfaces like soil and sand. Yet today we spend most of our time standing, walking and running on unnatural, hard surfaces.

Natural relief



Before (left) After (right)

Clinical research supports Vasyli Medical's success in restoring the foot's natural (neutral) alignment, enabling the feet, ankles and legs to function dynamically as they were intended. The unique orthotic is built into our entire collection of footwear.

Orthotics - A Natural Cure...

It is estimated that the majority of the population over pronates due to the hard, flat unnatural surfaces we walk on daily.

Over pronation often disrupts normal knee function and hip alignment and increases forces on the lower back muscles.

These biomechanical deficiencies are common in all types of patients from children to the elderly, from top athletes to people with a sedentary lifestyle.

By realigning the lower limb to its natural angle Vasyli® orthotics ensure correct foot function and help relieve common biomechanical complaints.

In combination with the practitioner's regular treatment regime, orthotics are usually successful in alleviating chronic pains.



the Problem...



the Cure

Ideal for:

- Plantar Fasciitis
- Metatarsalgia
- Tired/Aching Legs
- Lower Back Pain
- Ilio-tibial Band Syndrome
- Achilles Tendonitis
- Hallux Abducto Valgus
- Patello Femoral Pain
- Tibial Stress Syndrome

“ Supported by science, proven by customer satisfaction. ”

The Vasyli Think Tank™ is at the very heart of Vasyli.

Vasyli®+Dananberg 1st Ray Orthotic

Treating 1st Ray Function Just Got Easier!


Developed with podiatrist Howard Dananberg, renowned as a world leader in concepts relating to 1st ray function, the Vasyli Dananberg orthotic is a fast and simple alternative to a 1st Ray cutout. The 2 removable plugs allow you to customize the orthotic to assist in maximizing 1st ray function. In addition, this orthotic is also heat moldable to personalize the rearfoot control. Sold as pairs. Latex free.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
5630-11-01	XS	3½ - 5	4½ - 6
5630-11-02	S	5½ - 7	6½ - 8
5630-11-03	M	7½ - 9	8½ - 10
5630-11-04	L	9½ - 11	10½ - 12
5630-11-05	XL	11½ - 13	—
NEW! 5630-11-06	XXL	13½ - 15	—

NEW Size XXL coming in June 2014



Sorbon shock dot designed to dissipate force along the angle of strike.



Heat-moldable EVA allowing customization of the orthotic.

2 removable plugs allowing you to customize the orthotic further to assist in maximizing first ray function.



Howard J. Dananberg
Podiatrist & member
of the Vasyli Think Tank

Vasyli®+McConnell Slimfit High Heel Orthotic

Fitting orthotics into women's high heels and sandals is now much easier

Developed with physical therapist Jenny McConnell, renowned as a world leader in treatment of the lower limb and taping techniques. With a new slimmer profile and enhanced metatarsal support, the new Vasyli McConnell slimfit device "disappears" under the foot making it the perfect solution for hard-to-fit footwear.

- Enhanced metatarsal control
- Slim design fits most heeled shoes & sandals
- Dual-density compound
- Supports & realigns the foot

Item No.	Size	Women's Shoe Size
5629-89-01	XS	4½ - 6
5629-89-02	S	6½ - 8
5629-89-03	M	8½ - 10
5629-89-04	L	10½ - 12



Jenny McConnell
Physiotherapist & member
of the Vasyli Think Tank

The Vasyli Think Tank™ is at the very heart of Vasyli.



Brian Hoke

Received his degree in Physical Therapy from Indiana University and earned his post-professional Doctor of Physical Therapy from Boston University. He is board certified as a Sports Clinical Specialist by the American Board of Physical Therapy Specialties. Brian is co-owner and Director of Atlantic Physical Therapy, a physical therapist owned private practice in Virginia Beach, Virginia. Brian lectures extensively in the USA and internationally.



Vasyli®+Hoke Supination Control Orthotic

By utilizing a triplanar (three dimensional) perspective to normalize foot biomechanics, this orthotic effectively addresses forefoot valgus, plantarflexed 1st ray and forefoot equinus. The excessive supinator typically has:

- In the Sagittal Plane, a forefoot equinus frequently accompanied by a plantarflexed 1st ray
- In the Frontal Plane, a forefoot valgus
- In the Transverse Plane, a forefoot adductus with a prominent fifth metatarsal base.

The Vasyli Hoke device integrates triplanar features to address each of these anatomical considerations.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0816-09-957	XS	3½ – 5	4½ – 6
0816-09-965	S	5½ – 7	6½ – 8
0816-09-973	M	7½ – 9	8½ – 10
0816-09-981	L	9½ – 11	10½ – 12
0816-09-999	XL	11½ – 13	N/A



Thomas G. McPoil

Currently a Professor in the School of Physical Therapy at Regis University in Denver, CO. He has specialized in the evaluation and management of foot and ankle disorders since 1979. His scholarly contributions have systematically examined foot and ankle function from both scientific and clinical perspectives.



Vasyli®+McPoil Tissue Stress Relief Orthotic

During activity, a certain degree of foot pronation is necessary for normal lower extremity mobility. To accommodate for this individual specific degree of "normal" foot pronation and simultaneously prevent over pronation, the Vasyli McPoil Tissue Stress Relief Orthotic incorporates a dual density design that achieves both functional control and excellent compliance. The dual-density, total arch contact principle of the Vasyli McPoil foot orthotic design is applicable in the treatment regimes of both athletes and non-athletes with lower extremity conditions such as: anterior knee pain, medial tibial stress syndrome, posterior tibialis tendinopathy, chronic plantar heel pain ("plantar fasciitis/fasciosis"), general arch pain, and mid-foot pain.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0816-10-005	XS	3½ – 5	4½ – 6
0816-10-013	S	5½ – 7	6½ – 8
0816-10-021	M	7½ – 9	8½ – 10
0816-10-039	L	9½ – 11	10½ – 12
0816-10-047	XL	11½ – 13	N/A

Influencing product development to benefit both practitioner and patient alike.

Wider and thinner fore foot profile than the original Armstrong Orthotic!



Top view

Bottom view


Vasyli®+Armstrong II Sensitive Feet Orthotic

Whether you are dealing with an insensate foot, a hyper-sensitive foot type, or a diabetic foot, reducing pressure and shear forces should be playing a significant roll in your patients orthotic treatment regime.

Developed with podiatrist David Armstrong, renowned as a world leader in diabetic foot care, the Vasyli Armstrong II is the most advanced foot orthotic in its class. Featuring patented Glidesoft technology; the dorsal surface of the Vasyli-Armstrong II orthotic moves with the plantar surface of the foot dramatically reducing shear/frictional forces. This in-turn assists in the reduction of skin lesions and the potential for infection/ulceration.

The VAI is re-inventing the way leading health professionals care for sensitive feet! Ideal for the diabetic foot, patients suffering from rheumatoid arthritis or plantar fat pad atrophy causing painful sensitivity issues. This orthotic offers a clinically validated and scientifically tested shear and pressure reduction solution. Customize by using a standard heat gun and/or adding Vasyli postings. Sold as pairs. Latex free.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0815-80-406	XS	4-5	5-6
0815-80-414	S	6-7	7-8
0815-80-422	M	8-9	9-10
0815-80-430	L	10-11	11-12
0815-80-448	XL	12-13	13-14
0815-80-455	XXL	14-15	—



David Armstrong
Professor of Surgery
Director of Lower
Extremity Research and
the Southern Arizona Limb
Salvage Alliance (SALSA)
University of Arizona,
Tucson. Member of the
Vasyli Think Tank.



Top view

Bottom view

- The heat moldable high density shell (A) extends through the lateral aspect of the fore foot assisting in locking up the mid tarsal joint enhancing the propulsive phase of gait
- The PU fore foot (B) provides both comfort and function
- The 1st ray plantarflex zone (C) facilitates fore foot stability and Hallux dorsiflexion
- The PU shock dot in the heel (D) assists in attenuating heel strike forces

Vasyli®+Prior Sports Orthotic

The Vasyli + Prior incorporates a subtle 4mm heel raise—perfect for helping with tight posterior muscles and associated conditions that are often exacerbated in flat lasted sports footwear. The rear foot post of 3° + the 4mm raise forms a smart medial wedge for control without the “bulk.” Low profile, medial arch contour ensures comfort; the high density, heat moldable shell provides further customization, enhancing fit and performance.

The Vasyli Prior Sports Orthotic embodies the essential elements of function and control while satisfying the fit characteristics required in low profile / hard to fit performance footwear. Basically any boot of shoe with a stud, cleat or spike.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0815-80-463	XS	3½ – 5	4½ – 6
0815-80-471	S	5½ – 7	6½ – 8
0815-80-489	M	7½ – 9	8½ – 10
0815-80-497	L	9½ – 11	10½ – 12
0815-80-505	XL	11½ – 13	—
NEW! 0816-29-419	XXL	13½-15	—

NEW Size XXL coming in June 2014



Trevor Prior
Director of Premier podiatry
Limited and P2L, Trevor currently
consults as a Podiatric
Surgeon at Homerton University
Hospital in the UK. Member
of the Vasyli Think Tank.

Designed for tailor-made biomechanical alignment.

Vasyli® Custom Easyfit Insoles

Ideal for multiple shoe styles. This model features a lateral cut-away providing an easy fit into hard-to-fit footwear, including men's moccasins, loafers, and ladies (low-heeled) court shoes. Sold as pairs. Latex free.



Vasyli® Shock Absorber Orthotic

Has a mild intrinsically posted rearfoot to help reduce excess pronation and is made from soft PU and shock absorbing materials to cushion the foot. Whilst not as corrective as the standard Vasyli custom orthotics it does offer mild control not normally found in shock absorbing insoles. Built-in metatarsal raise, supporting the forefoot area. Great for supinators and can be posted using any of the Vasyli postings and additions. Sold as pairs. Latex free.



Full Length Only

Vasyli® Custom Extended Slimfit Insoles

Specifically designed to fit ladies shoes, such as fashion shoes and high-heels. The flexible base allows easy contouring of the shank of the shoe. The rearfoot angle is reduced to 2 degrees due to the high heel factor (tibio talar axis) whose action inverts the foot as the heel is raised. A built-in metatarsal raise ensures correct alignment of the metatarsal bones, helping to prevent Metatarsalgia and Morton's Neuroma. Sold as pairs. Latex free.



Sulcus Length Only



Green = Low Density

Available in 3/4 Length & Full Length

Ideal for elderly, pediatric and arthritic patients who require special consideration. Patients adapt easily to the mild functional control while the correct biomechanical position is maintained. Sold as pairs. Latex free.



Blue = Medium Density

Available in 3/4 Length & Full Length

Designed for the "average" patient with general orthotic needs, blue density orthotics provide excellent biomechanical control. Sold as pairs. Latex free.



Red = High Density

Available in 3/4 Length & Full Length

Designed for high performance, this range provides both durability and solid biomechanical control. Red, high performance orthotics are ideal for athletes, high activity users and overweight patients. Sold as pairs. Latex free.

Green		Blue		Red							
3/4 Length	Full Length	3/4 Length	Full Length	3/4 Length	Full Length	Easyfit	Shock Absorbers	Extended Slimfit	Size	Mens Shoe Size	Womens Shoe Size
5609-34-01	5609-33-01	5609-29-01	5609-27-01	5609-30-01	5609-28-01	—	—	—	KS	Youth 12–13	Youth 12–13
5609-34-02	5609-33-02	5609-29-02	5609-27-02	5609-30-02	5609-28-02	5609-31-01	—	—	KL	1½ – 3	1½ – 4
5609-34-03	5609-33-03	5609-29-03	5609-27-03	5609-30-03	5609-28-03	5609-31-02	5629-94-05	5609-32-01	XS	3½ – 5	4½ – 6
5609-34-04	5609-33-04	5609-29-04	5609-27-04	5609-30-04	5609-28-04	5609-31-03	5629-94-01	5609-32-02	S	5½ – 7	6½ – 8
5609-34-05	5609-33-05	5609-29-05	5609-27-05	5609-30-05	5609-28-05	5609-31-04	5629-94-02	5609-32-03	M	7½ – 9	8½ – 10
5609-34-06	5609-33-06	5609-29-06	5609-27-06	5609-30-06	5609-28-06	5609-31-05	5629-94-03	5609-32-04	L	9½ – 11	10½ – 12
5609-34-07	5609-33-07	5609-29-07	5609-27-07	5609-30-07	5609-28-07	5609-31-06	5629-94-04	—	XL	11½ – 13	—
5609-34-08	5609-33-08	5609-29-08	5609-27-08	5609-30-08	5609-28-08	—	—	—	XXL	13½ – 15	—

Vasyli® Heel Cups

Specifically designed to provide instant relief from painful heel conditions. The Heel Piece is also designed for hard-to-fit footwear such as ladies fashion shoes and men's and women's dress shoes where other orthotics may be too bulky. Sold as pairs. Latex free.

5609-67-01 Small, Fits shoe size: Men's 5-8; Women's 6-9

5609-67-02 Large, Fits shoe size: Men's 9-13; Women's 10+



Vasyli® Orthotic Wedges and Additions

A full range of self-adhesive Vasyli orthotic wedges and additions have been designed to assist in creating a better position for the foot structure. The fact is all feet are not the same and often there's a need to modify the orthotic for specific problems that require extra biomechanical control. Package of 10. Latex free.


5609-39 Metatarsal Raise, pkg. of 10

Wedges			Degree
Forefoot		Rearfoot	
Valgus	Varus	Varus	
—	5609-37-01	5609-35-01	2
5609-36-01	5609-37-02	5609-35-02	4
5609-36-02	5609-37-03	—	6

Heel Lift			Lift
Small	Medium	Large	
5609-38-01	5609-38-04	5609-38-07	4 mm
5609-38-02	5609-38-05	5609-38-08	6 mm
5609-38-03	5609-38-06	5609-38-09	8 mm

Mixed Addition Pack. This set provides you with a sampling of the Vasyli additions. Includes 1 pair of heel lifts, forefoot varus, rearfoot varus, forefoot valgus and metatarsal raisers. † **5629-95**

†Size and degrees vary per pack.



Vasyli® Medical Starter Packs.

Purchase everything you need to get started making custom insoles, for your patient population, in one convenient kit. **All kits include a free heat gun and manual.**



Active Starter Pack.

Includes 4 pairs of Dananberg insoles (sizes xsmall through large); 3 pairs full-length blue insoles (sizes small, medium and large); 2 pairs full-length red insoles (sizes large and xlarge); 1 pkg. Rearfoot Varus 2 degree (item #56093501) and 2 packs of mixed additions pkg (item #562995).

0816-12-167 Active Starter Pack

All Rounder Starter Pack.

Includes 6 pairs of Dananberg insoles (2 pr of each sizes small through large); 3 pairs 3/4-length blue insoles (sizes small, medium and large); 2 pairs Extended SlimFit insoles (sizes xsmall and small); and 2 packs of mixed additions pkg (item #562995).

0816-12-175 All Rounder Starter Pack

The Sports Starter Pack.

Includes 5 pairs of Dananberg insoles (sizes xsmall through xlarge); 5 pairs of Prior insoles (sizes xsmall through xlarge) and 2 packs of mixed additions pkg (item #562995).

0816-12-183 Sports Starter Pack

Tools



Economy Magnetic

Inclinometer. Same tool that is used in Vasyli workshops is now available for your convenience.

0815-23-570



Vasyli Heat Gun.
0815-23-547

Vionic® Sandals for Men

The ideal motion control footwear for after-sport recovery.

Durable VIBRAM outsole for excellent traction on various surfaces. Biomechanically contoured, motion-control footbed made from triple-density EVA for stability in the rear foot and cushioning in the heel and forefoot. EVA midsole and VIBRAM outsole. Received APMA Seal of Acceptance.

VIONIC™
WITH
ORTHAHEEL®
TECHNOLOGY



Vionic® Bryce. ⚠

Performance men's toe post sandal with flexible leather uppers, padded jersey liner and a soft toe post for a comfortable fit. Available in whole sizes only, medium (D) width. Latex free.

Black/Grey	Khaki	Men's Shoe Size
0815-89-324	0815-89-399	7
0815-89-332	0815-89-407	8
0815-89-340	0815-89-415	9
0815-89-357	0815-89-423	10
0815-89-365	0815-89-431	11
0815-89-373	0815-89-449	12
0815-89-381	0815-89-456	13



Vionic® Boyes. ⚠

Performance men's adjustable back-strap sandal offers secure comfort for active lifestyles. Flexible leather uppers feature fully adjustable straps with hook-and-loop closure and soft padded jersey liner. Available in whole sizes only, medium (D) width. Latex free

Black	Taupe	Men's Shoe Size
0815-96-535	0815-96-592	8
0815-96-543	0815-96-600	9
0815-96-550	0815-96-618	10
0815-96-568	0815-96-626	11
0815-96-576	0815-96-634	12
0815-96-584	0815-96-642	13



Vionic® Base. ⚠

Provides the foundation for stable sports recovery. The slip-on, toe post sandal is built upon a triple-dense midsole, and features a pop of color etched in architectural detailing. Available in whole sizes only, medium (D) width. Latex free.

Black/Silver	Choc/Red	Men's Shoe Size
0815-96-329	0815-96-386	8
0815-96-337	0815-96-394	9
0815-96-345	0815-96-402	10
0815-96-352	0815-96-410	11
0815-96-360	0815-96-428	12
0815-96-378	0815-96-436	13

⚠ This product contains natural rubber latex or dry natural rubber.



Vionic® Traverse. ⚠

This slide-on, adjustable sandal crosses conventional boundaries to provide the comfort, support and stability demanded by athletes of all levels. Available in whole sizes only, medium (D) width. Black/Silver only. Latex free.

Item No.	Men's Shoe Size
0815-96-444	8
0815-96-451	9
0815-96-469	10
0815-96-477	11
0815-96-485	12
0815-96-493	13

Half sizes for men should go up a size (i.e. 10½ goes to 11)

Vionic® Sandals for Women

9

The ideal motion control footwear for after-sport recovery.

VIONIC™
WITH
ORTHAHEEL®
TECHNOLOGY



Vionic® Cascade. ⚠

Performance women's toe post sandal with soft leather uppers, padded jersey liner with bungee strap overlay and a soft nylon webbing toe post. Available in whole sizes only in medium/wide (B/C) width. Latex free.

Black/Blue	Chocolate/Orange	Khaki/Red	Grey/Yellow	Women's Shoe Size
0815-89-464	0815-89-530	0815-89-605	0815-89-670	5
0815-89-472	0815-89-548	0815-89-613	0815-89-688	6
0815-89-480	0815-89-555	0815-89-621	0815-89-696	7
0815-89-498	0815-89-563	0815-89-639	0815-89-704	8
0815-89-506	0815-89-571	0815-89-647	0815-89-712	9
0815-89-514	0815-89-589	0815-89-654	0815-89-720	10
0815-89-522	0815-89-597	0815-89-662	0815-89-738	11



Vionic® Mojave. ⚠

Performance women's sandal with flexible leather uppers and two-tone webbing incorporates an elasticized toe ring plus hook-and-loop adjustable strap for a customized fit. Available in whole sizes only in medium/wide (B/C) width. Latex free.

Charcoal/Berry	Grey/Yellow	Women's Shoe Size
0815-96-162	0816-15-079	6
0815-96-170	0816-15-087	7
0815-96-188	0816-15-095	8
0815-96-204	0816-15-103	9
0815-96-212	0816-15-111	10



Vionic® Muir. ⚠

Performance women's back-strap sandal with flexible leather uppers and two-tone webbing incorporates hook-and-loop adjustable strap in both the forefoot and midfoot. The design offers a customized fit. Available in whole sizes only in medium/wide (B/C) width. Latex free.

Charcoal/Blue	Taupe	Women's Shoe Size
0815-96-220	0815-96-279	6
0815-96-238	0815-96-287	7
0815-96-246	0815-96-295	8
0815-96-253	0815-96-303	9
0815-96-261	0815-96-311	10



Vionic® Atlas. ⚠

Athletic style and support have evolved into the Atlas, the ideal toe post sandal for recovery after sports or a long day on foot. Features built-in orthotic footbed and dual density midsole for support and comfort. Available in whole sizes only in medium (B) width. Latex free.

Black/Silver	Grey/Coral	Women's Shoe Size
0815-95-966	0815-96-014	6
0815-95-974	0815-96-022	7
0815-95-982	0815-96-030	8
0815-95-990	0815-96-048	9
0815-96-006	0815-96-055	10

Half sizes for women should go down a size (i.e. 7½ goes to 7)



Vionic® Sprint. ⚠

Sprint to recovery with ease in this slip-on, after-sports sandal. With a pop of color for shelf appeal and femininity, plus medial side gore for individualized comfort and a triple-dense midsole, it rewards the foot with support and stability at every finish line. Available in whole sizes only in medium (B) width. Latex free.

Black/Silver	White/Citron	Women's Shoe Size
0815-96-063	0815-96-113	6
0815-96-071	0815-96-121	7
0815-96-089	0815-96-139	8
0815-96-097	0815-96-147	9
0815-96-105	0815-96-154	10

⚠ This product contains natural rubber latex or dry natural rubber.



VIONIC™
WITH
ORTHAHEEL®
TECHNOLOGY



Vionic® Tide II. ⚠

Fun and comfortable, the favorite Tide flip flop sandal makes it easy to support your feet and add color to your wardrobe as the weather warms up. Available in Women's whole sizes in three fashionable colors! Latex free.

Black	Black Sequin	Bronze Metallic	Bronze Sequin	Pewter Metallic	Gold Sequin	Women's Shoe Size
0815-96-659	0816-260-27	0815-96-725	0816-260-92	0815-96-790	0816-261-67	5
0815-96-667	0816-260-35	0815-96-733	0816-261-00	0815-96-808	0816-261-75	6
0815-96-675	0816-260-43	0815-96-741	0816-261-18	0815-96-816	0816-261-83	7
0815-96-683	0816-260-50	0815-96-758	0816-261-26	0815-96-824	0816-261-91	8
0815-96-691	0816-260-68	0815-96-766	0816-261-34	0815-96-832	0816-262-09	9
0815-96-709	0816-260-76	0815-96-774	0816-261-42	0815-96-840	0816-262-17	10
0815-96-717	0816-260-84	0815-96-782	0816-261-59	0815-96-857	0816-262-25	11



Tide Sequins coming in Spring 2014

Vionic® Tide Sequins Sandal. ⚠

The icon of Orthaheel® technology is now offered with a sequined upper and a subtle dot design on the footbed for traction. Geometric patterned tread on the outsole. Women's medium/wide width, whole sizes only. Latex free.



Vionic® Ryder. ⚠

A classic flip flop sandal with built-in orthotic support and added adjustability. Ryder also features a padded jersey-lined upper and a soft nylon webbing toe post for even greater comfort. Available in Men's whole sizes in two classic colors. Latex free.

Black/Grey	Chocolate/Tan	Men's Shoe Size
0815-96-865	0815-96-931	7
0815-96-873	0815-96-949	8
0815-96-881	0815-96-956	9
0815-96-899	0815-96-964	10
0815-96-907	0815-96-972	11
0815-96-915	0815-96-980	12
0815-96-923	0815-96-998	13



Vionic® Unisex Wave Flip Flops. ⚠

The popular, iconic flip flop sandal from Vasyli features the built-in biomechanical Orthaheel footbed. The lightweight "heel" supports your patients' feet while helping to realign their lower legs and improve posture. Durable rubber outsole with wave-patterned tread helps improve traction. Latex free.



Unisex Wave Flip Flop			Men's Shoe Size	Women's Shoe Size
Black	Chocolate	Khaki		
5647-09	0815-53-668	5646-99	5	6
5647-10	0815-53-676	5647-00	6	7
5647-11	0815-53-684	5647-01	7	8
5647-12	0815-53-692	5647-02	8	9
5647-13	0815-53-700	5647-03	9	10
5647-14	0815-53-718	5647-04	10	11
5647-15	0815-53-726	5647-05	11	12
5647-16	0815-53-734	5647-06	12	13
5647-17	0815-53-742	5647-07	13	14

⚠ This product contains natural rubber latex or dry natural rubber.

Half sizes for men should go up a size (i.e. 10½ goes to 11)
Half sizes for women should go down a size (i.e. 7½ goes to 7)

Contact your Sales Representative for new styles in spring and fall!

Built-in orthotic support provides relief from common foot, ankle and knee pain

VIONIC[™]
WITH
ORTHAHEEL[®]
TECHNOLOGY



Vionic® Relax Women's Slippers ▲.

Soft and comfy women's slipper with terry uppers and hook-and-loop adjustable straps. This style can be worn year-round for much-needed support around the house. Fabric-covered TPR outsole with wave-pattern tread designed for indoor use. Built into the slipper is our award-winning bio-mechanical orthotic footbed which supports your feet and helps realign the lower legs—improving posture and balance. EcoFresh® antibacterial, odor-reducing footbed and linings. Latex free.

Black	Tan	Zebra	Women's Shoe Size
0815-80-588	0815-48-833	0816-17-455	6
0815-80-596	0815-48-841	0816-17-471	7
0815-80-604	0815-48-858	0816-17-489	8
0815-80-612	0815-48-866	0816-17-497	9
0815-80-620	0815-48-874	0816-17-505	10



Vionic® Olivia. This ballet flat dances circles around all others. The richly textured man-made uppers are complete with contrast lining and a pretty little bow across the toe. Olivia may look dainty, but she's perfectly constructed to provide support. Available in whole and half sizes (except sizes 10 and 11); medium (B) width only. Latex free.

Black Patent	Women's Shoe Size
0816-18-172	6
0816-18-180	6.5
0816-18-198	7
0816-18-206	7.5
0816-18-214	8

Black Patent	Women's Shoe Size
0816-18-222	8.5
0816-18-230	9
0816-18-248	9.5
0816-18-255	10
0816-18-263	11

Vionic® Rhythm.

This lightweight fitness walker with breathable mesh keeps you moving in comfort and style. Features our biomechanical orthotic technology, supportive midsole brace and durable rubber outsole with flex grooves. Available in whole and half sizes (except sizes 10 and 11); medium (B) width only. Latex free.



Black/Teal	White/Blue	Red/Grey	Women's Shoe Size
0816-18-370	0816-18-479	0816-28-825	6
0816-18-388	0816-18-487	0816-28-833	6.5
0816-18-396	0816-18-495	0816-28-841	7
0816-18-404	0816-18-503	0816-28-858	7.5
0816-18-412	0816-18-511	0816-28-866	8
0816-18-420	0816-18-529	0816-28-874	8.5
0816-18-438	0816-18-537	0816-28-882	9
0816-18-446	0816-18-545	0816-28-890	9.5
0816-18-453	0816-18-552	0816-28-916	10
0816-18-461	0816-18-560	0816-28-924	11

Half sizes for women should go down a size (i.e. 7½ goes to 7)

▲ This product contains natural rubber latex or dry natural rubber.

Contact your Sales Representative for new styles in spring and fall!




Orthaheel® Relief 3/4 Length.

EVA orthotic with Tri-Planar Motion Control Technology. Cushioned shock dot in the heel for added shock absorption. Cut-away on the outside edge making it easy to fit into most mens and womens flat or low heeled shoes, or shoes without removable insole. EcoFresh® anti-bacterial odor absorbing cover. Contoured around the heel and arch areas to achieve 100% foot contact.


Orthaheel® Relief Full Length.

EVA orthotic with Tri-Planar Motion Control Technology. Cushioned shock dot in the heel for added shock absorption. Can be trimmed in forefoot if necessary. 4-degree rear foot wedge to provide support and control which helps prevent over-pronation. EcoFresh® anti-bacterial odor absorbing cover. Contoured around the heel and arch areas to achieve 100% foot contact.


Orthaheel® Active Full Length.

EVA orthotic with re-enforced hardened plastic. Cushioned shock dot in the heel for added shock absorption. Can be trimmed in forefoot if necessary. Designed to provide greater control in faster paced activities such as running and fast walking. 4-degree rear foot wedge to provide support and control which helps prevent over-pronation. EcoFresh® anti-bacterial odor absorbing cover. Contoured around the heel and arch areas to achieve 100% foot contact.

Relief 3/4 Length	Relief Full Length	Active Full Length	Size	Men's Shoe Size	Womens Shoe Size
0815-93-300	0815-93-250	0815-93-201	XS	3.5-5	4.5-6
0815-93-318	0815-93-268	0815-93-219	S	5.5-7	6.5-8
0815-93-326	0815-93-276	0815-93-227	M	7.5-9	8.5-10
0815-93-334	0815-93-284	0815-93-235	L	9.5-11	10.5-12
0815-93-342	0815-93-292	0815-93-243	XL	11.5-13	N/A

Orthaheel® Gel Heel Device.

Made from soft PU Gel. Designed to fit almost all types of footwear. Great shock-absorbing qualities with EcoFresh® ant-bacterial deodorizer top cover. Softer rear-foot wedge provides low-level control of over-pronation

Gel Heel Device	Size	Men's Shoe Size	Women's Shoe Size
0815-93-359	Small	5-7.5	6-8.5
0815-93-367	Medium	8-10.5	9-11.5
0815-93-375	Large	11-13	12-14.5



2014

LOWER LIMB BIOMECHANICS WORKSHOPS



BROUGHT TO YOU BY

PATTERSON MEDICAL & VASYLI MEDICAL

In keeping with our educational based philosophy, Vasyli Medical has the pleasure in announcing the schedule for the 2014 biomechanics and orthotics workshops.

This year our Biomechanical workshops are being presented by Brian Hoke, Howard Dananberg and Tom McPoil. All are highly regarded in the lower limb biomechanics field and have many years of practical experience to share with you.

If your practice sees patients with common lower limb conditions such as plantar fasciitis, shin splints, patella femoral syndrome, and low back pain then these are the perfect workshops for you. Our courses are ideal if you have been considering dispensing orthotics and you are perhaps a little apprehensive as to how to "get started" and the appropriate techniques to use when evaluating for orthotics. On the other hand, if you already use orthotic devices, are looking for a refresher on the biomechanics and relevant theories, or want to learn more about refining your prescriptions then these straight forward approaches will be of great interest to you.


These workshops are a great blend of theory and include plenty of practical "hands on" attendee participation. Most importantly the content is geared towards making better patient outcomes a reality for your facility.

We look forward to seeing you there!

Please check www.vasylimedical.com/seminars for updates



Vasyli Medical Orthotic Patient Education Brochure
Package of 25
0815-28-348



Vionic Sandal Flyer
Package of 25
0816-14-148



Vionic Sandal Poster
(22" x 28")
0816-14-122



Vionic 4-Fold Family Product Flyer
Package of 25
0816-15-145









Orthaheel Insole Flyer
Package of 25
0816-02-770







Grid Display
Versatile, can be used as a floor stand or a wall mount display.
*0816-12-563

* Item number reflects our 2013 Spring Kiosk. See your Patterson Medical sales representative for new footwear selections.



Vasyli Medical Product Matrix

CATEGORIES	Armstrong Sensitive Feet	Dananberg 1st Ray Function		Hoke Supination	McConnell High Heels	McPoil Tissue Stress	Prior Sports	
								
Red Box indicates application								
CONDITION								
Back Pain								
Knee Pain								
Foot/Heel Pain				Arch/Forefoot				
Diabetic								
USAGE								
Running								
Walking								
Day to Day								
Lateral Activity (Basketball etc.)								
After Activity/Early Morning								
SHOE TYPE								
Athletic						Low Profile		
Men's Casual	Depth Essential							
Women's Casual	Depth Essential							
Men/Ladies Casual (Narrow)								
Women's Dress (High Heel)								
Men's Dress								
Custom Fit Shoes								
Kids Shoes								
Football Boots		Deep Heel In Shoe	Deep Heel In Shoe				Low Profile	
Work Boots								
Work Boots (High Heat Environment)								
SIZES AVAILABLE	XS, S, M, L, XL	XS, S, M, L, XL,	XXL	XS,S,M,L,XL	XS,S,M,L,XL	XS,S,M,L,XL	XS,S, M,L,XL	XXL
PRODUCT AVAILABLE	NOW	NOW	JUNE	NOW	NOW	NOW	NOW	JUNE

CATEGORIES	Custom Range 3/4 & Full Length	Easy Fit	Extended Slimfit	Shock Absorber	Heel Pain Stabilizer	Vionic® Footware
						
Red Box indicates application						
CONDITION						
Back Pain						
Knee Pain						
Foot/Heel Pain						
Diabetic						
USAGE						
Running						
Walking						
Day to Day						
Lateral Activity (Basketball etc.)						
After Activity/Early Morning						
SHOE TYPE						
Athletic						
Men's Casual						
Women's Casual						
Men/Ladies Casual (Narrow)						
Women's Dress (High Heel)						
Men's Dress						
Custom Fit Shoes						
Kids Shoes						
Football Boots			Shallow Heel		Shallow Heel	
Work Boots	Full Length					
Work Boots (High Heat Environment)						
SIZES AVAILABLE	KS,KL,XS,S,M,L,XL,XXL	KL,XS,S,M,L,XL	XS,S,M,L	XS,S,M,L,XL	S,L	5-13 US
PRODUCT AVAILABLE	NOW	NOW	NOW	NOW	NOW	NOW



VIONIC™
WITH
ORTHAHEEL®
TECHNOLOGY

Partnership Performance Innovation

Our feet were designed to walk on natural, soft surfaces like soil and sand. When the feet strike soil or sand, these surfaces give way to accommodate the lower limb's natural 6 degree angle. Unfortunately, we give our feet unnatural hard, flat surfaces like pavements and floors. This forces our feet to roll over and our arches to flatten – a condition called excess pronation. Poor alignment of the feet and legs can cause wear and tear to other parts of the body, often disrupting normal knee function and hip alignment as well as increasing forces on muscles in the lower back.

Vionic Group products help restore natural foot function.



Learn more from the leaders in biomechanics and foot orthotics with Vasyli Medical's innovative mobile application.



Available to Rehab Clinicians exclusively through Patterson Medical


PATTERSON
MEDICAL
www.pattersonmedical.com

Call toll free:
1-800-323-5547
In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
In Canada fax 1-877-566-8611

©2014 Patterson Medical Products, Inc.
All rights reserved. Printed in USA.
081587476LIT

VIONIC®