

Recovery Series

ENCOMPASS POWERTOWER™ CLINICAL COMPLETE PACKAGE

Distributed by
PERFORMANCE HEALTH

Formerly Patterson Medical
7101459

CC Package accessories: SCRUNCH,
Weight Bar, Leg Pulley System and Press Bar

Also includes standard unit accessories:
CCP and BAPS Adapter, Retractable Dip Bars,
Squat Handle Bars, Slide Distance Regulator



PRODUCT SPECIFICATIONS

RESISTANCE	1% to 72% of bodyweight
INCLINE LEVELS	26 calibrated levels, with infinite adjustment between levels
IN USE	108" x 38" x 80" (L/W/H) [2.7 m x 1.0 m x 2.0 m]
FOLDED	20" x 38" x 65" (L/W/H) (.5 m x 1.0 m x 1.7 m)
UNIT WEIGHT	218 lbs (99 kg)
STORAGE	Hydraulic lift assists folding Rolls upright
USER HEIGHT	Up to 6'10" (2.1 m)
CAPACITY	Maximum total 650 lbs (295 kg)
POWER	110~220 VAC; 50~60 Hz
CONSTRUCTION	Steel reinforced extruded aluminum rails
WARRANTY	Frame - 5 years; Parts & Upholstery - 1 year; Foam & Rubber - 90 days; Motor - 2 years

The Encompass PowerTower™ allows incremental load changes during exercise, providing the utmost versatility for post-injury or post-surgery rehabilitation as well as sports specific and athletic training. Motorized level adjustment improves client access to glideboard and ease of transitions between exercises.

The Clinical Complete Package adds four additional accessories to provide comprehensive support to clinicians in the utilization of the Encompass PowerTower – from rehabilitation to sport-specific power training and progress assessments.



SCRUNCH
For early closed-chain functional movement patterns in a partially loaded environment.

WEIGHT BAR
Provides attachment point for Biomechanical Ankle Platform System (BAPS board).

LEG PULLEY SYSTEM
Provides control and safety by enabling you to set the rolling distance of the glideboard.

PRESS BAR
Provides a solid handhold during squat exercises.