

Sit to Stand Over Ground On a Treadmill



PG-360 Features	
Weight Capacity	360 lbs.
Overall Height	7' 9"
Max Patient Height	7' 2"
Max Height on 6" TM	6' 8"
Max Yoke Height	7' 9"
Powered Yoke Adjustment	28"
Base Width	35.4" / XL 39"
Inside Base Width	30.7" / XL 34.3"
Leg Height	8"
Leg Length	48"
Handlebar Range	22.5"
Caster Diameter	5"
Stool	Quick Release PhysioStool
Harness Included	Size L 42.5" - 46.8"
Warranty	3 Yr Frame,
	3 Yr Actuator,
	1 Yr Harness,
	1 Yr Electronics,
	3 Month Battery
Optional Harness Sizes	XXS: 24.4" - 29.5"
	XS: 29.9" - 33.1"
	S: 33.5" - 36.2"
	M: 36.6" - 41.3"
	L: 42.5" - 46.8"
	XL: 47.2" - 51.2"



PhysioGait® Product Overview

Dynamic Unweighting Systems

1 Flex-Yoke

Flexible Yoke Adjustments

Mini-Squats & Variable Unweighting

3 Actuator

Powerful & Quiet

Transfer from Sit to Stand

360° Yoke + 2

Flip and Store

Rotate & Lock in **Any** Position


Gait Scale 4

Gait Displacement Scale

Measures uni & bi lateral support

Displays

- Left & Right
- Total
- Differential



5 Battery

Long-Lasting

+ Wall-Mounted Charger

7 L or XL Base

Door Ways & Treadmills

Designed to Work in Your Facility

Harness 6

Ideal Lifting Location

Comfortable and Adjustable

Handles 8

Numerous Hand Grips

Vertical & Horizontal Handles

9 Physio-Stool

Magnetic Stool

Ease of Use for Therapists

For More Information Visit

PhysioGait.com