



Nitro Aluminum Rollator, 10" Casters Height Adjustable, Removable Back Support, Seat and Lever Locks




FIGURE A



FIGURE B



INSET C



INSET D



INSET E



INSET F



RETAIL BOX



RTL10266WT

- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button (Figure A)
- Back support height easily adjusts with tool-free thumb screw (Figure B)
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage (Figure C)
- Handle is located on seat
- Plastic clip holds the Nitro closed when frame is folded (Figure D)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (Figure E)
- Handsome and removable zippered storage bag Unique bag attachment keeps bag securely in place when rollator is open or folded (Figure F)
- Lightweight, aluminum frame
- Comes in retail box

PRODUCTS	Description
RTL10266	Red, 1/cs
RTL10266WT	White, 1/cs
RTL10266BK	Black, 1/cs

SPECIFICATIONS	Description
Handle Height	33.5" - 38.25"
Unit Dimensions	23" (W) x 27.75" (L)
Seat Dimensions	18" (W) x 20.5" (H) x 10" (D)
Weight	17.5 lbs.
Carton Shipping Weight	21.5 lbs.
Weight Limit	300 lbs.
Warranty Limited Lifetime	(frame)
5 Year Limited Warranty	(brake and handle housing)

ACCESSORIES	Description
892	Cane Holder, 1/ea
10266-CH	Nitro Cup Holder, 1/ea
10266-4OX	Oxygen Tank Holder - Standard, 1/ea
10266HD-4OX	Oxygen Tank Holder - Heavy Duty, 1/ea

PARTS	Description
1026602	1/ea
1026603	1/ea

PARTS	Description
1026601-RN	1/ea
1026601-RNR	1/ea
1026601-G	1/ea
1026601-GR	1/ea
1026601-W	1/ea
1026601-WR	1/ea
1026632	1/ea
1026630N	1/ea
1026604-B	1/ea
1026604-G	1/ea
1026605-R	1/pr
1026605-W	1/pr
1026605-G	1/pr
1026638	1/ea
1026611	1/PR
1026610-G	1/ea
1026610-R	1/ea
1026610-W	1/ea
1026633-R	1/ea
1026633-BK	1/ea
1026633WT	
1026625	1/ea
1026624	1/ea
1026625-AN	1/EA
1026624-AN	1/EA
1026617	1/ea
1026617-AN	1/EA
1026615	1/ea
1026616	1/ea
1026641	1/ea
1026642	1/ea
1026612-L	1/ea
1026612-R	1/ea
1026609-FSG	1/ea
1026609-FSW	1/ea
1026626-R	1/ea
1026626-G	1/ea
1026626-W	1/ea
1026644	1/ea
1026607-R	1/ea
1026607-G	1/ea
1026606-W	1/ea
1026613	1/pr
1026619	1/pr
1026631	1/ea
1026618	1/pr
1026618-AN	1/PR
1026614	1/pr
1026614-AN	1/EA
1026620	1/pr
1026622	1/ea
1026621	1/ea
1026621-AN	1/EA
1026623	1/pr