

INTRODUCING

The Rolyan Splint Bath

A Clean, Smart Design



PATENT-PENDING


Smart Splint Bath Technology Inside

Our state-of-the-art splint baths are the first of their kind. Made from a thermoset composite material known for its durability, insulation, flammability resistance, and aesthetic appeal, the Rolyan smart splint bath stands apart from traditional metal splint baths commonly found on the market. From a touch screen interface that provides the clinician with complete control over the unit's parameters, to an easily accessible transparent hinged lid, to a convection heat exchanger that provides quick and uniform heat distribution, the Rolyan smart splint bath covers it all.

- + Auto on/off settings
- + Water evacuation via pump/valve system
- + Precise temperature reading (+/-2° F)
- + convenient side storage container for scissors




MADE IN USA. PATENT-PENDING.



Systematic Approach to Infection Prevention

99.99% Kill Rate Achieved



Unit sanitation and patient cross contamination is a common concern in hospital and clinical environments. Our proprietary cleaning system uses both UV sterilization cycling and Sanitation Mode, a high-heat effective sterilization of the bath. This dual cleaning system insures the water is 99.99% free of microbes. In addition, the automatic drainage pump makes water evacuation quick and easy for tub cleaning and water replacement. This systematic approach makes the Rolyan splint bath the most sanitary on the market, without exception**.


- + Convection heat exchanger cycles water over UV light every 2 minutes
- + Nightly high heat Sanitation Mode gives an extra level of sterilization
- + Reduces the risk of cross contamination

**Jans, A. (2019) Antibacterial Activity and Efficacy of Kismet Holdings Group's Test Device Using a Custom Suspension Time-Kill Procedure. Round Rock, TX: Microchem Laboratory.




Intuitive User Experience

Traditional splint bath controls can be imprecise and clumsy to use. Our full-color touch screen allows users to browse and select preloaded thermoplastic settings from the top manufacturers or create their own. With an easy to read timer and temperature always on display, the clinician can focus on the patient and splinting. A full menu of additional features allows for clinician customization.




Home screen with factory presets of thermoplastic sheet recommendations.



Automatic water evacuation with the push of a button.



Prepare the bath to your ideal settings with a simple touch.



Custom presets to fit clinician preferences and non-programmed thermoplastics.



For more information, contact your Performance Health Sales Rep.

Performance Health
28100 Torch Pkwy | Suite 700 | Warrenville, IL 60555 | USA
Copyright ©Kismet Holdings Group, LLC 2019
Form #500.600 Patent-Pending