



People.Health.Care.

ReadyWrap®

Put your best foot forward with our diverse foot coverage options.



Low-stretch material
Mimics short-stretch bandaging — the gold standard of treatment.

50% overlapping straps
Helps prevent gapping.¹

Various foot choices
Coverage options for all presentations.



People.Health.Care.

ReadyWrap®
Adjustable
compression
wrap



Pair with ReadyWrap® Calf:

ExoAnklet™
Circular knit anklet
with 23–32 mmHg
compression.

ReadyWrap Foot SL
Low-stretch wrap
with slim foot and
adjustable ankle.

ReadyWrap Foot CT NEW
Low-stretch wrap
with great contour
and containment.

ReadyWrap Foot
Low-stretch wrap
with adjustable foot
and ankle.

Off-the-shelf	Yes	Yes	Yes	Yes
Edema/lymphedema management	Mild	Mild to Moderate	Moderate to Severe	Moderate to Severe
Adjustability	None	Ankle	Ankle and dorsum of foot	Ankle and dorsum of foot
Fits in everyday footwear	Yes	Yes	Some*	Some*
Foot length options	Regular, Long	Regular, Long	Regular, Long	Regular
Left/Right Specific	No, Unilateral†	Yes	Yes	Yes
Max. MTP circumference (A^{1c}/A^c)	33 cm (A ^{1c})	36 cm (A ^c)	41 cm (A ^{1c}) / 31 cm (A ^c)	Not applicable
Max. Ankle circumference (B^c)	41.5 cm	38 cm	43 cm	38 cm
Max. Heel circumference (Y^c)	43 cm	Not applicable	Not applicable	48 cm
Available colors	BLACK	BEIGE BLACK	BEIGE BLACK	BEIGE BLACK
Sold as	Pair (2)	Individually (1)	Individually (1)	Individually (1)
Includes a ReadyWrap Liner	No	Yes	Yes	Yes

Sources

1. Ehmann, S., and Bock, K. (Jan. 2018), “Case report to demonstrate the need for selection criteria for optimal adjustable Velcro wrap prescription”, The Journal of Advanced Wound Care, North American Supplement, vol. 27, no. 1

* Footwear which adjusts across the dorsum may be needed to accommodate adjustable straps.

† ExoAnklet is a unilateral design that can be worn on either foot, or used as one to wash and one to wear.

