



HERE COMES THE SUN

red light therapy wrap

extra long

This Extra Long Wireless Red Light Therapy Wrap combines light wavelengths of both red light (660nm) and infrared light (850nm) for targeted cellular-healing, and potential pain and inflammation reduction. This versatile, wireless, extra long, and flexy therapy wrap contains 450 hard working LED light chips to provide convenient, comfy light therapy for a larger and longer surface area.

Use this wrap all over the place – feet, ankles, knees, hips, low back, neck, shoulder, elbow, wrist, and hands, get creative.

Safe, effective, science-backed, affordable light therapy. Yep, for real.

why red light therapy?

Red Light therapy is used by medical professionals and studied by scientists for its ability to improve how the cells in our body produce energy. The long list of potential benefits continue to amaze people - including pain relief, reduced inflammation, increased blood flow and improved collagen production - just to name a few. Recovering from an injury? Wanting to improve your performance? Looking for safe, effective, enjoyable and easy multi purpose healing wrap? We made it just for you.

how to use

Place your wrap (skin contact) anywhere you need some help! Secure with the Velcro, or don't. Up to you. Turn it on. Take a break. Ahhhh. Convenient built-in 20 minute timer. Adjust intensity for sun-like warmth. 20 minutes a day, 3-7 times/week.



features

- Lightweight, convenient, and comfortable
- Adjustable Light Intensity settings - 100% Light, 50% Light
- 150 medical grade, water resistant 0.5 W LED Lights
- 3 light chips per diode – 450 LED chips
- More light for more healing.
- Combined Red Light 660nm and Infrared Light 850nm wavelengths
- Bendy, Flexible, Durable water-resistant neoprene wrap with an adjustable strap
- Cordless mobile battery powered recovery tool. No wires!
- Rechargeable Lithium Polymer Battery - 3-4 uses/charge

what's included?

- Red Light Therapy wrap
- Adjustable Velcro strap
- USB-C charging cord and USB Adapter (5V-2A)
- User guide
- 2 year limited warranty


pricing

msrp: \$269



tech specs

LED Lights	<ul style="list-style-type: none">• 75W Wrap – 150 x 0.5W, medical grade, LED Lights• 3 LED light chips/light- 450 total LED light chips for a large, uniform light coverage
Light Wavelength	<ul style="list-style-type: none">• 660nm (Red):850nm (NIR - Near-Infrared) wavelengths in a 1:2 ratio• 2 NIR light chips, 1 RED light chip
Light Irradiance or Intensity	<ul style="list-style-type: none">• 25-30 mw/cm2 maximum intensity at 0"• 'Sun-Like' intensity• Continuous wave of light - no pulsing
Light Use & Adjustability	<ul style="list-style-type: none">• Automatic 20-minute timer• 2 levels of light brightness - 100% or 50%
Battery	<ul style="list-style-type: none">• Li-polymer battery, 5000mah, 5V• 3-4 uses per charge• 2-3 hour charging time
USB-C Charging Adapter	<ul style="list-style-type: none">• Rated input- output 5V-2A
Size & weight	<ul style="list-style-type: none">• 31" x 6.3", useful for multiple body regions• 0.3 lbs, light and convenient
Material	<ul style="list-style-type: none">• Flexible and Durable water-resistant neoprene pad with an adjustable strap
Warranty & Light Lifetime	<ul style="list-style-type: none">• 2 Year Limited Warranty• Light Lifetime 50,000 hours
Safety – EMF & Light Flicker	<ul style="list-style-type: none">• EMF – 0.0mG at 0"• Light Flicker - Less than 1%, easy on the eyes
Conformance	<ul style="list-style-type: none">• FCC, RoHS
FDA Cleared	<ul style="list-style-type: none">• Drug-free pain relief



yesterday's medicine, today's science

At Fringe, our mission is a simple one: to create products, content, and community centered around healing. In a world that moves fast, exploring the healing traditions from our past offers us slower-paced, longer lasting solutions to our modern day health challenges.