



HERE COMES THE SUN

# red light therapy wrap

elbow, hand, wrist

The Fringe Elbow, Wrist & Hand Red Light Therapy Wrap combines light wavelengths of red light (660nm) and near infrared light (850nm) and an optional vibration massage setting.

Wrap your elbows, wrist, and hands for targeted cellular-healing, as well as potential pain and inflammation reduction. This ergonomic, wireless, and versatile therapy wrap contains 210 hard working LED light chips and 4 optional vibration massage areas to provide convenient, comfy, warm, massaging light therapy.

Safe, effective, science-backed, affordable light therapy. The perfect warm, healing, massage-y recovery tool for your busy hands, wrists and elbows.

### why red light therapy?

Red Light therapy is used by medical professionals and studied by scientists for its ability to improve how the cells in our body produce energy. The long list of potential benefits continue to amaze people - including pain relief, reduced inflammation, increased blood flow and improved collagen production - just to name a few. Recovering from an injury? Wanting to improve your performance? Looking for safe, effective, enjoyable and easy multi purpose healing wrap? We made it just for you.

### how to use

Place your wrap (skin contact) anywhere you need some help! Secure with the Velcro, or don't. Up to you. Turn it on. Take a break. Ahhhh. Convenient built-in 20 minute timer. Adjust intensity for sun-like warmth. Add gentle vibration massage - Ooohh. 20 minutes a day, 3-7 times/week.



## features

- Lightweight, convenient, and comfortable
- Adjustable Light Intensity settings - 100% Light, 50% Light
- 90 medical grade, water resistant 0.5 W LED Lights
- 3 light chips per diode – 270 LED chips
- More light for more healing.
- Combined Red Light 660nm and Infrared Light 850nm wavelengths
- Bendy, Flexible, Durable water-resistant neoprene wrap with an adjustable strap
- Cordless mobile battery powered recovery tool. No wires!
- Rechargeable Lithium Polymer Battery - 3-4 uses/charge

## what's included?

- Red Light Therapy wrap
- Adjustable Velcro strap
- USB-C charging cord and USB Adapter (5V-2A)
- User guide
- 2 year limited warranty


## pricing

msrp: \$229



## tech specs

<b>LED Lights</b>	<ul style="list-style-type: none"><li>• 45W Wrap – 90 x 0.5W, medical grade, LED Lights</li><li>• 3 LED light chips/light- 360 total LED light chips for a large, uniform light coverage</li></ul>
<b>Light Wavelength</b>	<ul style="list-style-type: none"><li>• 660nm (Red):850nm (NIR - Near-Infrared) wavelengths in a 1:2 ratio</li><li>• 2 NIR light chips, 1 RED light chip</li></ul>
<b>Light Irradiance or Intensity</b>	<ul style="list-style-type: none"><li>• 25-30 mw/cm2 maximum intensity at 0"</li><li>• 'Sun-Like' intensity</li><li>• Continuous wave of light - no pulsing</li></ul>
<b>Light Use &amp; Adjustability</b>	<ul style="list-style-type: none"><li>• Automatic 20-minute timer</li><li>• 2 levels of light brightness - 100% or 50%</li></ul>
<b>Vibration Massage Use</b>	<ul style="list-style-type: none"><li>• Optional vibration massage function with 100% light</li><li>• Approx. 150Hz at 2.3V.</li></ul>
<b>Battery</b>	<ul style="list-style-type: none"><li>• Li-polymer battery, 5000mah, 5V</li><li>• 3-4 uses per charge</li><li>• 2-3 hour charging time</li></ul>
<b>USB-C Charging Adapter</b>	<ul style="list-style-type: none"><li>• Rated input- output 5V-2A</li></ul>
<b>Size &amp; weight</b>	<ul style="list-style-type: none"><li>• 20" x 13", elbow &amp; wrist/hand</li><li>• 0.4 lbs, light and convenient</li></ul>
<b>Material</b>	<ul style="list-style-type: none"><li>• Flexible and Durable water-resistant neoprene pad with an adjustable strap</li></ul>
<b>Warranty &amp; Light Lifetime</b>	<ul style="list-style-type: none"><li>• 2 Year Limited Warranty</li><li>• Light Lifetime 50,000 hours</li></ul>
<b>Safety – EMF &amp; Light Flicker</b>	<ul style="list-style-type: none"><li>• EMF – 0.0mG at 0"</li><li>• Light Flicker - Less than 1%, easy on the eyes</li></ul>
<b>Conformance</b>	<ul style="list-style-type: none"><li>• FCC, RoHS</li></ul>
<b>FDA Cleared</b>	<ul style="list-style-type: none"><li>• Drug-free pain relief</li></ul>



yesterday's medicine, today's science

At Fringe, our mission is a simple one: to create products, content, and community centered around healing. In a world that moves fast, exploring the healing traditions from our past offers us slower-paced, longer lasting solutions to our modern day health challenges.