



ROLYAN ALIGNRITE WRIST SUPPORT WITH STRAP

The Rolyan AlignRite Wrist Support with Strap offers excellent support without limiting the user's finger or thumb movement. This support is constructed of woven, cotton-coated rubber fabric that features an open weave for ventilation. Adjust the notched stays to apply appropriate support and pressure.

Key Features:

- Features improved stitching for increased durability, wrap around wrist strap for extra support, and vertical two-piece design for an outstanding fit
- Sizing fits wrist circumference measures 4 ½" to 8 ¾" in size
- Constructed of latex-free woven, cotton-coated rubber fabric that is breathable and comfortable
- Wrap-around wrist strap offers additional support and features hook and loop strapping for a secure fit



SIZE	Wrist Circ.
XS	4½" to 5½"
Small	5½" to 6½"
Med	6½" to 7¼"
Large	7¼" to 8"
XL	8" to 8¼"



Scan to view
product online

SKU	Item Name	Specifications
081285048	Short - Left - X-Small	4.5 - 5.5"
081285063	Short - Left - Small	5.5 - 6.5"
081285089	Short - Left - Medium	6.5 - 7.25"
081285105	Short - Left - Large	7.25 - 8"
081285121	Short - Left - X-Large	8 - 8.25"
081285055	Short - Right - X-Small	4.5 - 5.5"
081285071	Short - Right - Small	5.5 - 6.5"
081285097	Short - Right - Medium	6.5 - 7.25"
081285113	Short - Right - Large	7.25 - 8"
081285139	Short - Right - X-Large	8 - 8.25"