



ROLYAN D-RING WRIST AND THUMB SPICA SPLINT ORIGINAL DESIGN - BEIGE

Rolyan D-Ring Wrist and Thumb Spica Splint allows healing after injury and helps decrease pain and inflammation of painful wrist and thumb tendonitis. Perfect for skier's thumb, gamekeeper's thumb and de Quervain's syndrome, and many other conditions that require rigid support such as ligament instabilities, sprains or muscle strains.

Key Features:

- Includes an aluminum stay that is pre-bent to 15 degrees of extension
- Non-elastic circumferential design evenly distributes pressure
- The D-ring closure system provides simple application and removal
- Sizing based on wrist circumference



Scan to view product online

SKU	Item Name	Specifications
081284223	Left, X-Small	Up to 5.75"
081284439	Right, X-Small	Up to 5.75"
081284678	Left, Small	5.75"-6.5"
081284777	Right, Small	5.75"-6.5"
081284801	Left, Medium	6.75"-7.5"
081284819	Right, Medium	6.75"-7.5"
081284827	Left, Large	7.75"-8.5"
081284835	Right, Large	7.75"-8.5"
081284868	Left, X-Large	Over 8.75"
081284231	Right, X-Large	Over 8.75"