

TREATMENT FURNITURE



PORTABLE PARALLEL BARS


Designed for ultimate versatility and unmatched durability, these Portable Parallel Bars provide a reliable training and rehabilitation solution for users of all levels. They deliver exceptional stability and support during use. The foldable design ensures effortless transport and a compact footprint for convenient storage. Made in the USA and shipped fully assembled.

Key Features:

- Welded aluminum frame and walking surface
- Adjustable, aluminum handrails
- Anti-microbial and anti-bacterial powder coat
- Anti-slip platform

Specifications

- Load capacity: 600 lbs on the base
- 300 lbs per bar
- Transport mode: 12.25-17.5"W x 60"L x 39.5"H, 55° Turn radius
- Extended mode: 45.5"W x 60"L x 28.5-41"H
- Adjustable handrail: 22.75-26.5"W
- Aluminum handrails
- Unit weighs 78 lbs
- Dual lock mechanism for height adjustment
- (4) 4" Total locking medical grade caster system



SKU	Item Name	Specifications
7031997	Portable Parallel Bars	22.75-26.5"W