

Material

Rolyan® TailorSplint™

Temperature

$\frac{1}{16}$ " (1.6mm) 150° to 160° F (65° to 70°C) water for 30 seconds

$\frac{3}{32}$ " (2.4mm) 150° to 160° F (65° to 70°C) water for 30 to 45 seconds

$\frac{1}{8}$ " (3.2mm) 150° to 160° F (65° to 70°C) water for 1 minute



Stretch

Moderate resistance

Time

$\frac{1}{16}$ " (1.6mm) 1 minute of working time

$\frac{3}{32}$ " (2.4mm) 2 to 3 minutes of working time

$\frac{1}{8}$ " (3.2mm) 3 to 5 minutes of working time

**LATEX
FREE**